

MED-MAX 8000

Introducing a **NEW**, Versatile, Clean and Simple Bariatric Bed



Bed weighs less than 300 lbs



FEATURES

- Sleeping Surface is comprised of 4 aesthetically pleasing, smooth sided tubular “box” sections that adjust in seconds from 35” to 39”, 42”, and 48” widths.
- Clear span under bed provides unobstructed access for lifts, etc.
- Bed travels straight up and down avoiding caster movement during operation.
- Standard baseboard bumpers eliminate wall damage
- Because the low position is only 9” the bed can also be used as a normal conventional High/Low bed
- For ease of portability the bed has a unique shipping and storage mode: Bed is comprised of 4 separate sections that bolt together in a package 1 1/2’ x 3’ x 4’ high that ride on its own casters
- Assembly of bed frame takes approximately 15 minutes
- Delivery by 1 person is easy, as the sections (all under 100 lbs) can be handled separately
- Four total lock, Hospital grade, 5” casters are standard
- Easy removable head and foot boards incorporate a unique shape for pendant placement.

Med-Mizer’s New Bariatric Low Bed
Implementing Versatility: A bed that can be used for both Bariatric and conventional Long Term Care applications

Options

- Battery Backup
- 4 Standard Colored Head and Footboards
- Trapeze
- 4 Inch Bed Extender
- 2 Position Side Rail that complies with all FDA Entrapment Zones
- Assist Bar

Specifications

- | | |
|-------------------------------|--------------------------|
| • Deck Height (low position) | 9” |
| • Deck Height (high position) | 25” |
| • Sleeping Surface | 35”, 39”, 42”, 48” x 80” |
| • Maximum Safe Working Load | 800lbs. |

MED-MAX 8000™

Offering the Ultimate Solution for Bariatrics

Adjustable Width



Lifetime warranty on all welds, 10 years on the frame, and 3 years on electrical components.

Approximate 15 Minute Assembly Time

- The Med-Max 8000™ offers ease of assembly .
- Eight standardized bolts promote quick set up.
- Four section design offers additional convenience for handling

Portability & Stand Up Feature

- One Person Delivery
- Fold-Up Feature
- Box-Tubular Frame Construction