


RETRACTABED[®] SYSTEM



Shown with embedded fixed control option
and four swivel casters

RetractaBed[®] SS Model

- Retracting Bed Frame - reduces potential falls and increases available area in room
- Embedded Controls - fixed controls in assist bar or head rail reduces replacement handsets and increases staff productivity
- Expanded Width - bed can expand to 42 inches with no tools
- Functional Lockouts - standard feature to prevent unnecessary operation
- Soft non-marking bumpers come standard to protect wall from damage
- The bed can be vertically transported and for ease of storage or getting through doors
- Clinical-Contour[®] - exclusive feature that reduces skin shear and the repositioning of residents



THE INNOVATIVE HEALTHCARE SOLUTION
www.med-mizer.com | 877-867-7365
Made in the U.S.A.

RetractaBed® SS Model Specifications

PHYSICAL SPECIFICATIONS

SS

- Deck Height Low 9"
- Deck Height High 25"
- Sleeping Surface 35" w x 80" L
- Length Retracted 71"
- Length Extended 85"
- Maximum Safe Working Capacity 500 lbs.
- Maximum Back Angle 64°
- Maximum Knee Angle 24°
- Mobility Height any height

ENVIRONMENTAL SPECIFICATIONS

- Operating Temperature 41° F to 104° F (5° C to 40° C)
- Operating Relative Humidity 20% to 90% @ 86° F / 30° C
- Operating Atmospheric Pressure 700 to 1600 hPa
- Storage Temperature -10° C to -50° C
- Storage Relative Humidity 20% to 90% @ 86° F / 30° C
- Storage Atmospheric Pressure 700 to 1060 hPa

ELECTRICAL SPECIFICATIONS

- Main Voltage Input 120 V AC
- Max Current Input 2.70 A
- Main Frequency Input 60 Hz
- Controller Output Voltage 24 V DC



Fall Risk
Non-Retractable Bed


Fall Prevention
RetractaBed®

LIMITED WARRANTY

- 15 year limited warranty on bed frame
- 5 years on drive system
- 1 year on wood products and casters
- Lifetime warranty on all welds

Approvals: IEC60601-1:2005 2nd edition, ANSI / AAMI ES60601-1:2005 2nd edition, and CAN / CSA-22.2 No 60601-1:2008

Endboard Options



Standard Board Shapes



ENGLISH OAK



MAHOAGANY



CHERRY



LIGHT OAK

*Custom boards are available



4" Bed Extender



¼ Length Side Rail



Floor Bumper Guard



¼ Length Side Rail with Fixed Controls



42" Bed Width Kit



Pivoting Assist Bars



Pivoting Assist Bars with Fixed Handset



Trapeze Assembly



Transport Dolly



THE INNOVATIVE HEALTHCARE SOLUTION

www.med-mizer.com | 877-867-7365

Made in the U.S.A.