

Falls PREVENTION Series



↓
100%

AllCare Bed and SafetyMat helped decrease bedfalls for an individual by 100% at St. Andrews Skilled Nursing Facility

EVEN PLANE

Create an even plane from the bed to the safety mat

SAFETYMAT AllCare MODEL

- **Bi-Fold Safety Mat** – 10" high safety mat creates an even plane from the top of the mattress to the safety mat. Turn a bed fall (impact) into a bed roll.
- **Reflective Handles** – added to enhance the safety of the caregiver and resident.
- **Advanced Falls Notification Package** – weight based sensor pad (reacts to weight) can work seamlessly with existing nurse call or independently with pager system and wireless door light.
- **Anti Skid Bottom** – added to bottom of the safetymat to reduce sliding.
- **Vertical or Horizontal Storage** – during "non" periods of rest store the safetymat in a convenient location in either a vertical or horizontal state.

Falls PREVENTION
Series

AllCare SafetyMat SPECIFICATIONS


SafetyMat Height..... 10"
 SafetyMat Width..... 30"
 SafetyMat Length..... 76"



Even Plane intended to reduce bedfalls by creating a wider width and turning an "impact" into a "roll"



"Light Reflective" handles to enhance caregiver safety



STORAGE (Vertical or Horizontal)



AntiSkid Bottom

Advanced Falls Notification System



Falls Notification Pad



Pager System



Wireless Door Light

PROUDLY MADE IN THE U.S.A.

(877) 867
MED-MIZER.COM



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS PREVENTION




BARIATRIC CARE



MOBILITY



REHABILITATION



WOUND CARE