

# FlexCHAIR™ Series

---



## FlexTilt® WC19 Transport Chair

- **WC-19 Certified** — meets standards to be utilized as a transport chair for residents up to 350lbs for various appointments
- **VBack Design®** — patent pending exclusive feature that keeps the resident centered in the chair while adding comfort and support while providing facilities the ability to standardize on a chair to accommodate more resident sizes
- **Arm Rest** — height adjustable arm rest provide flexibility to customize height as needed and as an added benefit the pivoting motion allows an easier transfer
- **Pressure Relief Strapping** — provide additional support and comfort by "flexing" to the individuals unique contour
- **Dual Steering Casters** — dual steering casters allow the staff to push the chair in a straight direction reducing back strain
- **Adjust-ability** — customize the chair to the right height and seat depth based on the residents needs
- **Made in the USA** — an affordable quality product with a commitment to keeping inventory and parts on hand.(no customs or long lead times)

## FLEXTILT CHAIR SPECIFICATIONS

<b>Frame</b> .....	Gauge steel 11ga - 16ga
<b>Degree of Tilt</b> .....	25°
<b>Degree of Recline</b> .....	30°
<b>Seat Width</b> .....	16", 20", 24"
<b>Seat Depth</b> .....	16.5"-18"
<b>Back Seat Height</b> .....	32"
<b>Seat Height</b> .....	19.5"-21.5"
<b>Weight Capacity</b> .....	350 LBS.
<b>Warranty</b> .....	3 year mechanical frame 1 year on components



### FlexChair Options



#### FlexTilt®

Available in 16", 20", 24"



#### FreedomFlex®

Available in 16", 20", 24"

PROUDLY MADE IN THE U.S.A.

(877) 867-7365

[MED-MIZER.COM](http://MED-MIZER.COM)



is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS  
PREVENTION




BARIATRIC  
CARE



MOBILITY



REHABILITATION



WOUND CARE